
Acupuncture Guide Manuel d'acupuncture

GB Research has shown that stimulating pens, such as Pain®Gone can also be used for non-invasive acupuncture. Please follow the simple guide to find your condition and locate the relevant points of acupuncture. When treating a condition for which more than three points are listed, start by trying out the first three points. If the effect is not satisfactory, vary the points until you find the best combination for you. Click 20 – 25 times on each acupoint you are treating. When treating migraine or headache, you must apply stimulation to the opposite side of the body. If the headache or migraine is on the right side of the head, you should apply stimulation to the left side of your body, for example on the left arm or the left hand.

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
Medi-Direct International Ltd
Unit 24, Wilford Industrial & Business Park
Ruddington Lane, Nottingham, NG11 7EP, UK
www.paingone.com

© Copyright, 2015, Medi-Direct International Ltd. All rights reserved. Copying or other reproduction of this document is prohibited without the prior written consent of Medi-Direct International Ltd.

Clinically tested
Testé cliniquement

Version No: 7 - 02/15