

OPERATING INSTRUCTIONS – ENGLISH

BEFORE USING THE STIMULATOR

- 1. INTRODUCTION 26
- 2. MEDICAL BACKGROUND 27
- 3. PRECAUTIONARY MEASURES..... 28

PATIENT INFORMATION

- 4. OVERVIEW 29
 - CONTROL BUTTONS 29
 - DISPLAY SYMBOLS 30
- 5. OPERATION 31
 - STEP-BY-STEP USE..... 31
 - AUTO STIMULATION MODE..... 33
 - TIMER 33
 - PROGRAM PAUSE..... 34
- 6. REPLACEMENT OF BATTERIES..... 35

PROFESSIONAL INFORMATION

- 7. PROGRAMS 36
 - PRESET PROGRAMS 36
 - CUSTOM PROGRAMS 38
 - PROGRAM LOCK 40
 - DISABLE THE AUTO FUNCTION..... 40
 - COMPLIANCE..... 40

ADDITIONAL INFORMATION

- 8. ACCESSORIES..... 41
- 9. CARE INSTRUCTIONS 41
- 10. TROUBLESHOOTING..... 42
- 11. FREQUENTLY ASKED QUESTIONS 43
- 12. TECHNICAL DATA 44
 - KEY TO THE SYMBOLS..... 45
- 13. ELECTRODE PLACEMENT GUIDE SEE COVER



1. INTRODUCTION

CEFAR Medical AB has been a medical supplier of professional TENS Pain Relief products for more than 30 years. The company was established in Sweden in 1975 and is currently represented in large parts of the world.

The Cefar stimulators are widely used and endorsed by medical professionals and top athletes around the world.


TENS is effective, has no side effects and it's economical. Through clinical research, areas of application for TENS and EMS are rapidly expanding. Cefar is working actively to further develop the method towards a natural treatment alternative for the consumer.

More information about TENS (Transcutaneous Electrical Nerve Stimulation), EMS (Electrical Muscle Stimulation) and our products can be found on our web site: www.cefar.se

CEFAR PRIMO PRO is a dual channel nerve stimulator (TENS) designed for use by both professionals and home users. It provides relief of a wide range of pains, like back pain and neuropathic pain. The stimulator features 7 preset pain programs and 3 custom programs. Program and amplitude can be individually set for each channel, which means that two different forms of stimulation may be used simultaneously. With CEFAR PRIMO PRO all kinds of TENS treatment are covered including mixed frequency stimulation.



4. OVERVIEW/CONTROL BUTTONS



1. ON/OFF

- Turns the stimulator on and off.
- Turns the stimulator off even when the key lock is activated.
- Can be used for terminating the stimulation at all times.

2. INCREASE (left and right channel)

- Increases the amplitude (intensity of stimulation). Press and hold the button to increase the amplitude continuously.
- **Note!** Always increase the amplitude cautiously.
- Terminates AUTO stimulation mode.
- Increases the number of minutes when setting the timer.
- Used when turning the program lock on and off (left button).
- Scrolls through choices in programming mode.

3. DECREASE (left and right channel)

- Decreases the amplitude (intensity of stimulation). Press and hold the button to decrease the amplitude continuously.
- Deactivates the key lock.
- Terminates AUTO stimulation mode.
- Decreases the number of minutes when setting the timer.
- Used when turning the program lock on and off (right button).
- Scrolls through choices in programming mode.

4. AUTO (left and right channel)

- Starts the AUTO test mode.
- Confirms the selected amplitude level in AUTO test mode.

5. PROGRAM (left and right channel)

- Selects a program (1-10).
- Pauses an ongoing program.
- Used when activating/deactivating the program lock (left button).

6. PROGRAMMING/CONFIRMATION (S=SET) (left and right channel). SEE CHAPTER CUSTOM PROGRAMS.

- Turns the stimulator into programming mode for custom programs 8-10 when pressed for 2 seconds.
- Confirms settings in programming mode.
- Initiates the timer setting when pressed during active stimulation.

7. CABLE HOLDER



7. PROGRAMS/PRESET PROGRAMS

No.	Program name	Frequency	Pulse duration	Treatment time	Preset time
1	Conventional TENS	80 Hz	180 μ s	No limit (minimum 30 minutes).	-- minutes
2	Burst TENS	2 Hz	180 μ s	20-45 minutes, 3 times a day. Remember that Burst TENS can cause muscle soreness.	-- minutes
3	Modulated pulse duration stimulation	80 Hz	70-180 μ s	No limit (minimum 30 minutes).	-- minutes
4	Treating sensitive areas/ Facial treatment	80 Hz	60 μ s	No limit (minimum 30 minutes).	-- minutes
5	Mixed-frequency TENS	80 Hz/ 2 Hz	180 μ s	Up to 3 hours per day.	-- minutes
6	Nausea treatment	10 Hz	180 μ s	About 10 minutes per session. Sessions can be given as often as every other hour or as little as three times a day.	10 minutes
7	CEFAR Flow TENS	80 Hz	70-180 μ s	No limit (minimum 30 minutes)	-- minutes



Each program is described with the indications it is best suited to treat, how the electrodes should be placed, and how long the treatment should last. To each example of indication is at least one Electrode Placement number, which refers to the suggested electrode placements shown in the cover of the manual.

Program description	Examples of indications/ Electrodeplacement no.
<p>Conventional TENS (high-frequency stimulation) is the first choice for both acute and long-term pain, both neurogenic and nociceptive. Conventional TENS is based on the Gate Control theory, which states that electric stimulation of A-beta fibers inhibits impulse transfer in the pain pathways.</p> <p>As a rule, the electrodes should be placed on or near the painful area, or over an area segmentally related to the painful area. Adjust the amplitude so that the stimulation gives strong, but pleasant paraesthesia - tingling. When using high-frequency stimulation, it's important to ensure that the patient has normal sensory of touch in the area where the electrodes are placed.</p>	<p>Neck pain3 Shoulder pain4 Elbow pain8 Rheumatic pain 10 Fractured rib 16 Lumbago 13, 14 Menstrual pain 17 Phantom limb pain 18 Hip pain 20 Osteoarthritic pain in the knee.. 22, 23 Wound healing..... 24</p>
<p>Burst TENS (low-frequency stimulation) is usually most effective for radiating (projected) pain in the arms and legs (rhizopathy), for conditions with reduced or changed sensory of touch, for deep muscular pain, or when the post-treatment effect of Conventional TENS is too short. Burst TENS treatment alleviates pain by stimulating muscles to release the body's own morphine-like substances, endorphins.</p> <p>Place the electrodes on a muscle in the painful area so that a visible contraction occurs, or on acupuncture points in the painful area. The stimulation should feel pleasant and give visible muscle contractions. Remember that the patient often feels the stimulation clearly before contractions become visible.</p>	<p>Mononeuropathy..... 11 Central pain 6, 7 Cervical rhizopathy.....5 Sciatica 19 Knee pain 23</p>
<p>Modulated pulse duration stimulation is a type of high-frequency stimulation where the pulse duration varies continuously. This can cause an undulating sensation, which may be more pleasant than a constant pulse duration. Use program 3 for pain alleviation and a massage effect on muscles like the trapezius.</p>	<p>Trapezius pain 12, 25 Lumbago 13, 14 For more examples, see Program 1.</p>
<p>Program 4 has a shorter pulse duration than the other programs. A short pulse duration is suitable for treating sensitive, nerve-rich areas like the face and high up on the neck. With a shorter pulse duration the amplitude can be increased, which makes it easier to find the right level without the patient experiencing pain.</p>	<p>Tension type headache1 Facial pain2 Neck pain3 Postherpetic neuralgia 15</p>
<p>Mixed frequency stimulation is also called Han stimulation (3 sec at 2 Hz and 3 sec at 80 Hz). Stimulation frequencies switch every three seconds, giving a combination of high- and low-frequency stimulation, which can offer a more effective treatment.</p> <p>Place the electrodes as you would for low-frequency stimulation - on a muscle in the painful area.</p>	<p>Neck pain3 Shoulder pain4 Lumbago 13, 14 Hip pain 20 Thigh pain..... 21</p>
<p>Program 6 is specially designed for treating nausea. Place the black electrode on acupuncture point PC6 (inside of the wrist) and the red one on acupuncture point TE5 (outside of the wrist).</p>	<p>Postoperative nausea.....9 Chemotherapy-induced nausea.....9</p>
<p>Always use 2 channels and 4 electrodes with Program 7. A modulated pulse duration stimulation is obtained, but since the channels are alternately active, a massaging, pumping effect is experienced.</p> <p>Place the electrodes as for high-frequency stimulation when a pain relieving and massaging effect is desired.</p>	<p>Trapezius pain 12, 25 Lumbago 14 Thigh pain..... 26 Shoulder pain4</p>



Programming
CEFAR PRIMO PRO
P8, P9, P10

Conventional (C)	Burst (B)	Pulse width/duration modulation (PWM)	Frequency modulation (FM)	Mixed frequency (MF)
<p>Frequency Preset: 80 Hz Range: 2-120 Hz</p>	<p>Frequency Preset: 2 Hz Range: 1-5 Hz</p>	<p>Frequency Preset: 80 Hz Range: 2-120 Hz</p>	<p>Pulse duration Preset: 180 µs Range: 50-300 µs</p>	<p>Frequency 1 (Burst) Preset: 2 Hz Range: 1-5 Hz</p>
<p>Pulse duration Preset: 180 µs Range: 50-300 µs</p>	<p>Time Preset: - - min. Range: 1-60, - - min.</p>	<p>Pulse duration 1 Preset: 70 µs Range: 50-300 µs</p>	<p>Frequency 1 Preset: 20Hz Range: 2-120 Hz</p>	<p>Frequency 2 Preset: 80Hz Range: 6-120 Hz</p>
<p>Time Preset: - - min. Range: 1-60, - - min.</p>		<p>Pulse duration 2 Preset: 180 µs Range: 50-300 µs</p>	<p>Frequency 2 Preset: 40Hz Range: 2-120 Hz</p>	<p>Stimulation time Preset: 3 s Range: 1-10 s</p>
	<p>Time Preset: - - min. Range: 1-60, - - min.</p>	<p>Modulation time Preset: 3 s Range: 1-10 s</p>	<p>Modulation time Preset: 5 s Range: 5-15 s</p>	<p>Time Preset: - - min. Range: 1-60, - - min.</p>
		<p>Time Preset: - - min. Range: 1-60, - - min.</p>	<p>Time Preset: - - min. Range: 1-60, - - min.</p>	



12. TECHNICAL DATA

CEFAR PRIMO PRO is a nerve stimulator (TENS) with two independent channels. It has seven preset programs and three custom programs. Programs and amplitude can be individually set for each channel.

Treatment with electrical stimulation requires the stimulation current to penetrate the resistance of the skin and the electrode, about 1500 ohms. CEFAR PRIMO PRO can penetrate this resistance and maintain a current of up to 60 mA. With a change in load from 100 to 1500 ohms, the stimulation current changes less than 10% from the set value.

The stimulator operates on either two non-rechargeable 1.5 V AA batteries or two rechargeable 1.2 V AA batteries, recharged in a separate charger.

CEFAR PRIMO PRO

Number of channels	2 (independent)
Constant current	Up to a resistance of 1500 ohm
Stimulation current/channel	0–60 mA
Waveform	Asymmetrical biphasic pulse, 100% compensated
Number of preset programs	7
Number of custom programs	3
Stimulation forms	Conventional Burst Modulated pulse duration/frequency Mixed frequency Alternated modulated pulse duration (CEFAR Flow TENS)
Max pulse duration	300 µs
Max frequency	120 Hz
Timer	0 to 60 min
Environment for storage, use and shipping	Temperature 10°C–40°C Air humidity 30%–75% Air pressure 700 hPa–1060 hPa
Power source	2 x 1.5 V AA non-rechargeable or 2 x 1.2 V AA rechargeable batteries
Current consumption for one channel, 80 Hz, 30 mA	150 mA
I r.m.s. max/channel	11 mA
Size	120 x 50 x 30 mm
Weight	approx. 180 g



KEY TO THE SYMBOLS



Read the operating instructions before to use.



Patient part type – Body Floating.



0413

Complies with MDD 93/42/EEC requirements. Notified body Intertek ETL Semko (0413).

INFORMATION RELATED TO ELECTROMAGNETIC COMPATIBILITY (EMC).

CEFAR PRIMO PRO is designed to be used in typical domestic or clinical environments and approved according to the EMC safety standard of EN 60601-1-2.

CEFAR PRIMO PRO emits very low levels in the radio frequency (RF) interval. Therefore it is not likely to cause any interference in nearby electronic equipment (radios, computers, telephones etc.).

CEFAR PRIMO PRO is designed to withstand foreseeable disturbances originating from electrostatic discharges, mains supply magnetic fields and radio frequency transmitters (such as mobile telephones).



